

Popliteal Positioner



Simplify Access for Popliteal Nerve Block

- Positioner for Easier Popliteal Access
- Foam with Antimicrobial cover
- Easy to Clean
- Low Slip Vinyl base

www.thinkmedicalinc.com

(919) 782-2485 (800) 849-8446 fax (919) 781-3064
thinkmedical@earthlink.net



Popliteal Positioner #1060 Instructions for Use: Place lower extremity on positioner to allow for optimal access within safe limits.

After each use clean thoroughly with Cavi-Wipes. Sani-Wipes/ similar hospital approved cleaner.